

SERRATUS™ VERTEBRAL BIOPSY NEEDLE*

NEW!

A cost-effective,
reliable, specialty
needle designed to capture a
quality vertebral core sample

* Designed for use with the
NeedleTech Bonus™ Access Needle

Features include:

Ergonomically engineered, radiolucent handles
Unique, proprietary Serratus™ cutting cannula
Blunt obturator (included with each needle)
Luer Lock attachment
FDA 510K approved



Exclusive Serratus™
Cutting Cannula

Specifications and offerings:

13 gauge and 15 gauge with 8 inch working length
Medical grade class VI approved materials
Offered packaged and sterile or bulk, non-sterile
Custom sizes and kits available upon request



NeedleTech Products, Inc. 81 West St. Attleboro, MA 02703
ph: 508-431-4000 www.needletech.com